



(MEA and HSSA tested for IVF laboratory use)

Features:

CO2 incubators, laminar flow hoods, glass, plastic, metals, equipment and ultra-sound probes.

Non – toxic, odourless, colourless and VOC free.

Rapid action against a wide range of Bacteria, Fungi & Encapsulated viruses

4th Generation Quaternary Ammonium Compounds.

Detailed Product Description	
Mode of Action:	Denatures membrane protein and phospholipid membrane.
Storage Temperature:	15 – 30 °C
Package:	500 ml Spray bottle
Shelf Life:	12 Months

For more information

Email: info@pylon-healthcare.com

Website: www.Pylon-Healthcare.com

Phone: +20 122 555 6925, +20 1000 922 377 , +20 38520295

Effective against:

- Staphylococcus aureus.
- Pseudomonas aeruginosa .
- Escherichia coli.
- Enterococcus hirae.
- Mycobacterium phlei (TbB).
- Salmonella typhimurium.
- MRSA - Bacteria (Methicillin resistant Staphylococcus aureus).

2 - Fungi:

- Aspergillus niger.
- Candida albicans.
- Trichophyton mentagrophytes.

3 - Viruses:

- Hepatitis B.
- HIV (AIDS - Virus).
- Rotavirus 1 minute.
- Norovirus (Calicivirus).
- Influenza A - virus (H5N1 - Bird flu).
- Influenza A - virus (H1N1 - Swine flu).
- SARS & MERS Cov.
- Picornavirus (Poliovirus - group) .
- Adenovirus (DNA - virus group) .